

Version History

November 2018



MultiSpec® Pro II

Version 1.5.1247

- Support of new LS-LAS (Laser) light source component
- Timer triggered Sample measurements
- Windows 10 (v1803) support

Version 1.4.1189

- SensoLogic (new ADK 2.5.1)
- SensoLogic provides advanced selection mode of SL applications
- OLUP new Camo engine V10.5 (versions between 10.2 and 10.5 not supported!)
- User-Input step enhanced: multiple fields per window
- AIC activation with user-defined settings
- Agro edition supports Auto DC
- S-PACT PEAXACT fixed; requires application server > v4.0.1

Version 1.3.1116

- Bug fixes: S.PACT PEAXACT v4

Version 1.3.1112

- Additional product edition: Agro
- New method-steps: User-Input, S.PACT PEAXACT v4, GPS
- Results extended: averaging filter, piece-wise fit
- Storage method-step features extended (tags, sub-directories)
- Migrated to .NET 4.6.1, VC runtime 2012
- Added support for Windows 10 (32 / 64 bit)
- Bug fixes, improved UI, stability and performance

Version 1.2.1041

- Result formula parser Create-Fit function added
- Value-Selector method-step added
- Storage method-step features extended
- Bug fixes, improved UI, stability and performance

Version History

November 2018



Version 1.1.983

- Support for Raman spectroscopy
- Support for Multi-Channel editing and channel views
- BEK Laser security status panel
- Bug fixes, improved UI, stability and performance

Version 1.0.773

- First Release
- 3 Product editions: Basic, Standard, Process