

Version History

December 2019

Windows Driver for PD-PCIe01

Version 1.0.4.0

- Fixes / optimization of I²C
- Deinitialization of scan controller modified to disable cleaning scans as well
- ISP (In system programming) supports new flash types (required for newer boards)

- Supports Windows 7 SP1 and Windows 10 (32 / 64 bit)
Windows XP and Vista no longer supported
- December 2019: New packages added for Windows 10
- New device class: "tec5 Data Acquisition Devices"

Version 1.0.0.0

- First official release

- Sleep-/Hibernation-Modes are not fully supported and should be disabled for unattended long-term measurements