

# Version History

November 2016



## Lens Colour Analyser TFM-1<sup>®</sup>

### Version history

#### **Version 2.0.1 (Update, available on our [website](#))**

This setup offers an update of an **existing** TFM-1 installation V2.x to the current version V2.0.1. For a complete new installation additional files are needed, which you find on the provided CD. Please check the parameter settings (TFM-1.ini) after finalizing the update. For an update of an older installation V1.x please contact the tec5 service.

#### **Version 2.0.1.0 (Nov. 2016)**

- Update and bug fixes for UV protection
- Scalability of the UV protection checks
- Supports Windows 7 / 8.1 / 10
- Bug fixes

#### **Version 2.0.0.0 (Nov. 2014)**

- Update EU Standard DIN EN ISO 8980-3
- Extensions for Australia/New Zealand AGS/NZS 1067 and US standard ANSI Z80.3
- MS SQL Server 2008 R2 Express SP2 (optional)
- custom-specific scanner modes (see manual)
- asymmetric a/b tolerance range
- Data base format Version 6
- VC runtime 2012
- Supports Windows XP / 7 / 8.x